

---

---

# Clarity

## *Getting Started*

ENG


---

---

Code/Rev.: M006/26A - 15 February 2008

Phone: +420 - 251 013 400  
Fax: +420 - 251 013 401  
[clarity@dataapex.com](mailto:clarity@dataapex.com)  
[www.dataapex.com](http://www.dataapex.com)

© **DataApex Ltd. 2008**  
Podohradská 1  
155 00 Prague 5  
The Czech Republic

Clarity<sup>®</sup>, DataApex<sup>®</sup> and <sup>®</sup> are trademarks of DataApex Ltd.  
Microsoft<sup>®</sup> and Windows<sup>™</sup> are trademarks of Microsoft Corporation.  
*DataApex reserves the right to make changes to manuals without prior notice.*  
*Updated manuals can be downloaded from [www.dataapex.com](http://www.dataapex.com).*



1	Brief Description.....	6
1.1	Hardware and software requirements.....	6
2	Installation .....	8
2.1	Software Installation.....	8
2.1.1	Windows Vista instalation.....	9
2.2	Hardware Installation .....	9
2.2.1	Hardware key (dongle) installation.....	10
2.2.2	The INT7 Card .....	10
2.2.3	Installation of the U-PAD Converter.....	11
2.3	Device setup and wiring.....	11
2.3.1	Standard cable for Clarity station.....	12
2.3.2	Chromatograph.....	12
2.3.3	Autosampler .....	14
2.4	Clarity configuration.....	14
3	Qualification procedures .....	17
3.1	Installation Qualification - IQ.....	17
3.2	Operational Qualification - OQ.....	18
3.3	Performance Qualification - PQ.....	18
4	Program structure and control.....	19
4.1	Main Clarity window.....	19
4.2	Instrument window .....	20
4.3	Method Setup dialog.....	21
4.4	Chromatogram window.....	22
4.4.1	Multi-detector chromatograms .....	24
4.4.2	Summary results table.....	24
4.5	Calibration window.....	25
4.6	Single Analysis dialog .....	26
4.7	Sequence window .....	27
4.8	Program exit.....	28
4.9	Data update and consistency .....	28
5	First “Analysis” .....	29
5.1	Program start .....	29
5.2	Create new project.....	30
5.3	Method development.....	31
5.3.1	Measurement and evaluation parameters .....	32
5.3.2	Create and assign a calibration .....	33
5.4	Data Acquisition.....	33



5.4.1	Sample description .....	33
5.4.2	Signal monitoring .....	34
5.4.3	Analysis run .....	35
5.4.4	Analysis monitoring .....	35
5.4.5	Analysis stop .....	36
5.5	Chromatogram evaluation .....	37
5.5.1	Display chromatogram .....	37
5.5.2	Chromatogram processing .....	38
5.5.3	Chromatogram modifications .....	39
5.5.4	Saving of a template method .....	40
5.6	Calibration .....	40
5.6.1	Opening a calibration file .....	41
5.6.2	Link the calibration file with chromatogram or template method.....	41
5.6.3	Calibration curve use .....	42
5.7	Printing.....	43
6	Connecting Autosamplers .....	44
6.1.1	AS + GC set – Active Sequence .....	44
6.1.2	AS + LC set – Active Sequence .....	47
6.1.3	AS - Passive Sequence (GC or LC).....	48
6.1.4	AS – Active Sequence + AS Control + A/D card .....	49
6.1.5	AS – Active Sequence + AS Control + digital acquisition .....	50
7	Troubleshooting.....	51
7.1	Locate your problem: .....	51
7.2	Problems at station start .....	52
7.2.1	DEMO – missing HW Keylock.....	52
7.2.2	DEMO – Keylock test failed .....	53
7.2.3	DEMO - WRONG S/N .....	53
7.2.4	DEMO - Trial Expired .....	54
7.2.5	DEMO (in the window header) .....	54
7.3	Problems upon collection of data .....	54
7.3.1	Data Acquisition – non-functional .....	54
7.3.2	Data Acquisition - Simulated .....	55
7.3.3	Other Error Messages .....	56
7.4	Hardware key .....	57
7.4.1	Rockey USB dongle (re)installation .....	57



- 7.4.2   Rockey LPT dongle (re)installation ..... 58
- 7.4.3   Sentinel dongle installation ..... 58
- 7.5    System Files (systeminfo.txt file) ..... 59
- 7.6    Sleep Mode ..... 60
- 7.7    Switching Users in Windows OS ..... 60



# 1 Brief Description

**Clarity** chromatography station is an effective tool for the acquisition, processing and evaluation of data from any gas or liquid chromatograph with analog output and from selected chromatographs with digital output.

In the maximum configuration, it is possible to measure on up to four chromatographs simultaneously, from which, each may be equipped with up to 12 detectors.

The station is equipped with support tools for automatic co-operation with chromatographs and autosamplers.

**Clarity** supports fulfillment of the requirements of **FDA's directive 21 CFR Part 11**.

The **Clarity** station automatically processes all data acquired using **CSW** stations.

## 1.1 Hardware and software requirements

The **Clarity** station runs in any language version of the following systems:

**Check the compatibility of your workstation:**

	INT5	INT7	INT9	U-PAD	U-PAD2	Net-PAD	CB11	CB20
<b>MS Windows NT(98/Me*)</b>	✓	✓		**		✓	✓	✓
<b>MS Windows 2000</b>	✓	✓	✓	✓	✓	✓	✓	✓
<b>MS Windows XP/Vista</b>		✓	✓	✓	✓	✓		✓

\* Support for Windows 98/Me ended with Clarity version 2.4.1.

\*\* U-PAD is functional only on MS Windows 98/Me, not on Windows NT.

**Note:** Only 32-bit versions of **MS Windows XP** and **Vista** are supported by the **Clarity** station.



### PC Configurations:

	Minimal	Recommended
<b>MS Windows Vista</b>	the same like MS Vista system requirements	
<b>MS Windows 2000/XP</b>	PC Pentium III/700 MHz, 256 MB RAM	PC Pentium 4/2 GHz, 512 MB RAM
<b>MS Windows NT/(98/Me*)</b>	PC Pentium/200 MHz, 32 MB RAM	PC Pentium/450 MHz, 128 MB RAM
<b>Monitor</b>	Resolution 1024x768, 64K (16 bit High color)	Resolution 1280x1024, 64K (16 bit High color)
	CD ROM drive, free slot for A/D converter (ISA/PCI/USB), USB or LPT for a dongle	

\* *Support for Windows 98/Me ended with Clarity version 2.4.1.*

Verify that you have:

- A free slot (or port) for your converter:  
 INT7, INT9, CB20      -    PCI slot  
 U-PAD, U-PAD2        -    USB port  
 Net-PAD                -    LAN connection

**Note:** *A full size PCI 2.0, 32 bit PCI slot is required. Low profile or PCI Express slots cannot be used.*

- A free **USB** slot or parallel (**PRINTER**) port for the hardware key (depending on the type).

**Caution!**

*For Windows NT the printer port hardware key is necessary, since Win NT does not support USB.*

- A CD ROM drive for software installation

**Note:** *When utilizing discontinued hardware (e.g. INT5 or CB11), consult the separate manual for requirements and compatibility issues.*



## 2 Installation


Verify that the package is complete (packing list).

**Caution!** *Install the Clarity before inserting any devices (INT7, U-PAD, etc.).*

### 2.1 Software Installation

Installation of the program requires 50 to 210 MB of free hard disk space depending on the number of components you want to install.


**Caution!** *Ensure you have Administrator access rights in your system before you proceed with the installation. Clarity users must have read/write access to the Clarity folders (C:\CLARITY and all subfolders).*

- Insert the CD into the drive.
- If the installation software does not start up on its own, select the INSTALL.EXE file and run it.
- The software installation wizard will take you through the installation procedure, including the creation of a **Clarity** item in the **Start - Program** menu and a **Clarity** icon  on the Desktop.

**Note:** *If you are only evaluating the program, do not enter the **User Code** and press the **Skip** button.*

- During installation of the drivers following message may occur: "Installation did not pass Windows Logo testing". If this occurs, select "Continue Anyway".





### 2.1.1 Windows Vista installation

When installing the **Clarity** software, **Windows Vista** may create a large number of warning screens (depending on the security levels set). These should be ignored for **Clarity** to function correctly.

**Note:** *It is highly recommended not to install the **Clarity** software into the PROGRAM FILES directory. Also, the **User Account Control (UAC)** in Windows Vista should better be disabled when working in the **Clarity**.*

**Caution!** *If you have already installed an older version of **Clarity** software than 2.6, unplug the hardware key during installation of the **Clarity** 2.6 version.*

## 2.2 Hardware Installation

The following chapters describe the installation of a protective **Hardware key** (dongle) and the brief installation of the most common integration converters **INT7** and **U-PAD**.

**Caution!** *Plug in the USB devices (ROCKEY USB hard lock or U-PAD) only after the installation of the **Clarity** station.*

A detailed description of the hardware and its installation including troubleshooting is available in separate manuals.



## 2.2.1 Hardware key (dongle) installation

The hardware key needs to be installed and present in the PC when using the Clarity Chromatography Station.

The drivers will automatically be installed during **Clarity** software installation. If the installation did not proceed automatically (older versions of MS Windows), refer to the chapter **7.4** - Hardware key on pg. **57**.

### Caution!

*Install Clarity software from the CD ROM first and only then insert the key into the PC.*

When you have the Sentinel hardware key, refer to chapter **7.4.3**.

### Printer port hardware key (dongle)

Plug the hard lock in the printer port between the printer port and any potential printer.

## 2.2.2 The INT7 Card

In **Windows 2000/XP/Vista** the drivers will be installed automatically during the installation of the **Clarity** software. Install the software before attaching the device to PCI slot.

### Caution!

*Install the **Clarity** before inserting any devices (INT7, U-PAD, etc.).*

- Turn off the **PC**
- Insert the **INT7** card into the **PCI** slot
- Start up the **PC**

### a) MS Windows 2000/XP and later

In **Windows 2000/XP/Vista** the drivers were installed automatically during the installation of the **Clarity** software.

### b) Previous versions of MS Windows

In other systems then **Windows 2000/XP/Vista** follow the procedure described in the separate INT7 manual.

- The driver is now installed, proceed to setup and wiring.

### Note:

*When installing multiple **INT7** cards, a PC restart is recommended after the installation of each card,*




in order to prevent the reassignment of INT7 acquisition channels between Clarity **Instruments**.

### 2.2.3 Installation of the U-PAD Converter

**Note:** In **Windows 2000/XP/Vista** the drivers will be installed automatically during the installation of the **Clarity** software. Install the software before attaching the device to USB slot.

- Install **Clarity** software from the CDROM.
- Connect the **U-PAD** with a cable to a **USB** port on the computer.
- After connecting the converter, **Windows** will recognize the new **Plug and Play** device and will attempt to install a driver. The **Found New Hardware Wizard** will appear:
- Select: "No, not this time."




- Select: "Install the software automatically".
- The rest of the installation will be carried out automatically. In the lower left corner of the screen, in the **Start** menu, select **Settings – Control Panel**; selecting the **System** icon, verify that the **Device Manager** tab contains a correctly installed "**Universal Serial Bus Controllers**" - "**U-PAD v. 1.3**" item.

## 2.3 Device setup and wiring

The standard **Clarity** station package consisting **A/D converter** (**INT7**, **U-PAD**, **Net-PAD**..) includes a **cable** set composed of signal, starting and digital output for connecting the **Clarity** station to the chromatograph.



### 2.3.1 Standard cable for Clarity station



- **Signal cables**

Labelled “**DET1**” to “**DET4**” (according to the number of channels), the cables are supplied as standard without connectors with only stripped, tinned endings – red (+), white (-) and shielding (analogue ground).

- **Starting (marker) cables**

Labelled “**IN1**” to “**IN4**” (according to the number of channels), ended with CINCH connector. One cable with free leads [red (+), shielding (digital ground)] for connection directly to the chromatograph or valve and one cable ended with a button for cases where starting contact is not available and it is necessary to perform start manually are supplied for each starting cable.


- **Cables of the digital outputs**

Relay contacts labelled, “**OUT 1R**” to “**OUT 4R**” (according to the number of channels), are ended with free leads. They are used for synchronizing autosamplers in the active sequence without an **AS Control** module.

- At the **converter end**, the cable is always ended with CANNON SUB D 37-pin connector.

### 2.3.2 Chromatograph

Connect cables according to one of the following diagrams in **Fig. 1**. Use symmetrical connection only in the case that you are sure that the chromatograph/detector is equipped with symmetrical output – it is necessary to read through the instructions for the corresponding chromatograph.



**Fig. 1. Connection of Clarity with chromatograph**

All current DataApex A/D Converters **INT7**, **INT9**, **U-PAD**, **U-PAD2** and **Net-PAD** use the same standard **INT7 Connector**.

**Note:** The **INT5** cable (obsolete type) requires a **reduction** to the INT7 card, which can be ordered separately.

Detailed description of individual connectors wiring is provided in the **Reference Guide** in the chapter **Technical Specification**.

#### **Principles of connection of signal cables:**

Signal inputs of the A/D converters are symmetrical: + (red), - (white), analogue earth (copper braiding).

#### **Caution!**

*Shielding must be connected. It serves not only as the shielding, but also as the analogue earth against which measurement takes place. In the case of asymmetrical output of a detector (only two connectors) shielding must be connected to the white lead! No lead of the signal cable may remain unconnected.*

Try to connect to the output of the chromatograph detector with the largest possible level of signal, usually indicated as **INTEGRATOR** (signal approx. 1V). The level of the signal on the output marked **RECORDER** is usually only approx. 10mV.



For easier changing of the wiring we supply **SV8 Terminal board** with screw contacts for **INT7** and **U-PAD A/D Converters**.

#### **Connection of starting cables:**

Starting input reacts to a change of the TTL logical level (5V) or to a connection by any contact (button, contact of relay).

Input implicitly reacts to a change from **HIGH** to **LOW** (or closing of a contact). The input function may be changed by switching the **Down** item of **Ext. Start/Stop** section from the **Method Setup - Measurement** dialog accessible from the **Instrument** window using the **Method - Measurement** command.


### **2.3.3 Autosampler**


The most typical autosampler connections are described in chapter **6** on pg. **44**.

The autosamplers controlled directly using an **AS Control** module (p/n A26) are described in separate manuals.

## **2.4 Clarity configuration**

**Note:** *Clarity is automatically preconfigured, however it is recommended to review the following chapter to assign custom detector names and signal units.*

- Start the **Clarity** station with the  icon on the desktop.
- Invoke the **System Configuration** dialog using the **System – Configuration** command.



**Fig. 2. System Configuration dialog**

- Set the number of instruments in the **Number of Instruments** field ①.

**Note:** A larger number of instruments can be set than the amount you have purchased. You will not be able to measure on the surplus instruments (indicated by a blue symbol of the curve on the tab), but you may use them e.g. for “offline” preparation of methods.


- Set the corresponding type of chromatograph (GC/LC/EA/GPC...) for each instrument with the **Instrument Type** switch ②.

From the **Setup Control Modules** list of installed equipment on the left drag the equipment according to your configuration into the selected Clarity **Instrument x** tab in the right-hand section of the dialog using the mouse or the **→** button ③.

**Note:** If necessary, add further equipment to the list on the left using the **Add** button ④ in the lower left-hand section of the dialog.



- The configuration dialog of the corresponding card or equipment is displayed with a double-click on its icon or using the **Setup** button:



**Fig. 3.** *Dialog for setting the converter INT7*

- Check that the correct number of channels is displayed, and, if necessary, perform further settings (e.g. set the names of the detectors).

**Note:** *You can change signal units using this dialog. More accurate description can be found in the manual of your A/D converter.*

- Press the **OK** button.





## 3 Qualification procedures

The quality of analytical data is an issue that has been gaining increased attention in many laboratories these days. One of the requirements for ensuring the reliability of generated results is the validation of all instrumentation and procedures that are used to acquire data.

For chromatography data stations usually three levels of validation (qualification) are relevant:

### 3.1 Installation Qualification - IQ

**Installation Qualification (IQ)** is a procedure confirming that the software was successfully installed and that the installation contains all the files of the correct version. Installation qualification is an integral part of the **Clarity** software installation procedure.

#### How to use Installation Qualification?

- Install the **Clarity** station according to the instructions of the **Installation Wizard**.
- After the installation has been completed, open the Windows **Start** menu and go to the **Programs – Clarity – IQ Report** item.
- The **IQ** window will be opened.
- If the installation has been correctly performed the status should read: "*Installation Qualification Test: Passed*".





Fig. 4. IQ window



- If the **Installation Qualification** fails, it is recommended to uninstall and then re-install **Clarity**. If it fails again contact DataApex user support.

**Note:** *The most common reason for obtaining a "Failed" result of the qualification is due to the installation of the station upgrade over an existing full version of **Clarity**. This itself does not produce any errors but since some of the files are preserved from the original installation, the checksums will not fit.*

- The **Installation Qualification** report can then be printed, copied to the MS Windows Clipboard or sent as an email.

## 3.2 Operational Qualification - OQ



**Operational Qualification (OQ)** is a procedure confirming that the data station is performing according to the manufacturer's specification. **The Operation Qualification is provided by Validation kit**, which consists of a precise peak generator and a set of methods and reports used in the validation process. **SST module**, an optional extension of **Clarity**, is also necessary.

The **Validation kit** (p/n: **CVK**) can be purchased separately, as well as **SST Extension** (p/n: **A22**).

## 3.3 Performance Qualification - PQ

**Performance Qualification (PQ)** is a procedure confirming that the analytical system is fit for a given type of analysis. Usually the overall system performance is tested by this procedure with respect to the requirements of the desired application. The evaluation of the system performance in **Clarity** is addressed using the dedicated Extension **System Suitability Test (SST)**.

The **SST module** (p/n: **A22**) can be purchased separately.



## 4 Program structure and control


The **Clarity** software has a hierarchic structure.

After start-up the main **Clarity** window will be displayed with the symbols of the configured instruments. After clicking on the chromatograph picture and entering the **User Name** (more information on **User Names** can be found in the **Reference Guide**) the **Instrument** window will be displayed. This window is used for acquisition and data processing on the connected chromatograph.

**Note:** *The **Clarity** station works with so-called Instruments. All detectors connected to one Instrument share a common time base.*


### 4.1 Main Clarity window

The main **Clarity** window is designed to set the station's configuration, select access rights and basic directories for saving data.



**Fig. 5. Clarity window**

	From
<b>Detectors</b>	
AA Detector 1	Int7 A/D Card driver
<b>GCs</b>	
GCs	HP5890 GC Driver
<b>AS</b>	
Sampler 1	HP7673 Sampler

The **System Configuration** dialog (opened using the  icon or **System – Configuration** command) is where number of connected chromatographs, their names, the selection of displayed symbols and the type of directly controlled equipment are all established.




By setting access rights (💡 or **System – User Accounts** command) you will allocate each user a name and password and you will determine the extent of authority for individual operation types (e.g. access to files, possible modification of calibration, integration parameters, etc.). Each user can set his/her own station appearance. For help on how to set a User Account, see Clarity Reference Manual.

## 4.2 Instrument window

The **Instrument** window is used for measuring and evaluating an analysis from a selected chromatograph. The window is displayed by clicking on the symbol of the relevant chromatograph in the station's main window. Depending on the number of the instruments, up to four independent **Instrument** windows can be displayed.

**Note:** *As many **Instruments** must be set and configured in the **Clarity** station, as you want to measure independent analyses.*



**Fig. 6. Instrument window**

Each **Instrument** window contains an information table, status line and analysis-processing



diagram. Windows are distinguished by line color in the analysis-processing diagram and Instrument name in the header.


Windows enabling performance of required actions on a given Instrument can be displayed by clicking on individual icons in the diagram or on commands from the menu.

### Caution!

*Certain Clarity windows (e.g. File Open or parameter setting dialogs) are modal. If a **modal window** is invoked you will not be able to work with any other Clarity window until you close it. It may seem as if the station does not react to your commands. Use the **MS Windows** shortcut **Alt** + **TAB** to switch to the modal window and close it.*

## 4.3 Method Setup dialog

Template method can be edited in the **Method Setup** dialog invoked from **Instrument** window using the corresponding command from the **Method** menu. The parameters determining and describing conditions for measurement and evaluation are set in the method and saved in the **template method file**.



**Fig. 7. Tabs of Method Setup dialog**

The **Method Setup** dialog is divided into the following sections accessible through icons directly from the **Instrument** window.

### Note:

*Some tabs (e.g. AS, GC, LC) may not be present depending on the configuration of the Instrument.*

**Event Table** – with events triggering variety of actions (digital outputs, run program, etc.)

**AS** – optional section for direct autosampler control



**Measurement** – section describing measurement conditions and containing optional setting of measurement duration



**Acquisition** – section of parameters for signal measurement – voltage range, sampling rate, ...



**Integration** – section with an integration table



**Calculation** – section with parameters for setting the type of calibration calculations






**LC/GC**– optional section for direct chromatograph control


**Advanced** – section with advanced parameters for chromatogram calculations

**Note:** *Individual instruments may add their own tabs into the method setup window. Have a look at the **List of Controlled Instruments** on the **DataApex** website ([www.dataapex.com](http://www.dataapex.com)). If you cannot find your instrument there, please let us know – the most frequently requested instruments would be developed first.*

## 4.4 Chromatogram window




The central window for display, modification and evaluation of chromatograms. This window can be displayed by clicking the  icon. Use  icon to display one of the stored chromatograms.




**Fig. 8. Chromatogram window**

The **Result Table** for a selected chromatogram can be found on **Results** tab. If it is not displayed automatically, click on the label of the **Results** tab in the bottom-left corner of the window or use command **Results - Result Table**. Modifications and evaluation of a chromatogram can be found in Chapter 5.5.

### Caution!




The **Chromatogram** window (and other **Clarity** windows) uses **Splitters** dividing the window into panes (graph, table, etc.). A double-click around the splitter (the cursor changes to ) enlarges the graph across the whole window. To restore both panes drag the **Splitter** to the desired location with mouse or use the **View - Show Both** command.

**TIP:** While working with enlarged part of chromatogram, small orientation graph may be useful. It can be displayed using **Display - Preview graph** command.



#### 4.4.1 Multi-detector chromatograms

In the case of multi-detector measurement, all signals are saved in one chromatogram.



Individual signals from detectors are displayed as independent curves. To work with individual signals, it is necessary to select the required signal from the **Chromatogram** menu or by using the coloured symbol in the **Overlay** toolbar. Each signal is independently worked with, just as with a single-detector chromatogram in the **Overlay** mode.

#### Caution!

Remember that the **Measurement** and **Calculation** tabs in the **Instrument** window, and also the selection of a calibration file, are common for all detectors. On the other hand, **Results**, **Performance** and **Integration** tabs belong to each detector independently.

#### 4.4.2 Summary results table


You will make use of **Summary Table** especially for reviewing and processing of results from sequences. Here, you can easily display all results from several chromatograms at once.


Summary Table											
	Sample ID	Injection Volume	1,2 cis DCE			1,2 dichloroethan			Benzen		
			Reten. Time [min]	Area [mV.s]	Response	Reten. Time [min]	Area [mV.s]	Response	Reten. Time [min]	Area [mV.s]	Response
	SR1_0112	0,500	6,624	603,357	603,357	7,620	662,980	662,980	8,227	1620,724	1620,724
	SR1_0211	0,500	6,620	569,703	569,703	7,614	646,687	646,687	8,223	1623,888	1623,888
	SR1_0410	0,500	6,624	624,349	624,349	7,618	660,908	660,908	8,227	1672,367	1672,367
	SR1_0411	0,500	6,650	606,219	606,219	7,646	652,067	652,067	8,255	1637,654	1637,654
	SR1_0612	0,500	6,630	629,497	629,497	7,622	678,584	678,584	8,227	1689,639	1689,639
	SR1_0710	0,500	6,600	394,932	394,932	7,590	553,746	553,746	8,199	1083,685	1083,685
←											
Results	Summary	Integration	Measurement Conditions	SST Results							






**Note:** Special add on module **SST** (p/n: **A22**) is available for statistical processing of results and monitoring of selected parameters.




Multiple chromatograms can be opened simultaneously in the Overlay Mode. To switch it on/off use the **File – Overlay Mode** menu command or  icon.

## 4.5 Calibration window




The window for creating, modifying and displaying calibration curves is displayed by clicking on the  icon.

The **Calibration Summary Table** and possibly also the calibration standard chromatogram are displayed on the **Compounds** tab.



**Fig. 9. Calibration window**


By switching to the individual compounds tabs you will display their calibration curves.



**Fig. 10. Individual compound tab**

**Note:** The procedure for creating a simple calibration and its use is included in Chapter 5.6 on pg. 40.

## 4.6 Single Analysis dialog




The dialog for defining individual analysis contains fields for sample description, sample naming wizard and buttons for controlling the analysis.


The screenshot shows the 'Single Analysis - (cb\_20\_init)' dialog box. It has two main sections: 'Analysis' and 'Control'. The 'Analysis' section contains fields for Sample ID (std 3), Sample (calibration standard level 3), Amount (0), Dilution (1), ISTD Amount (0), and Inj. Volume [ml] (0). There is a checkbox for 'Calibration Standard' which is checked, and a 'Level' dropdown set to '3'. A 'Method...' button is also present. The 'Control' section contains buttons for 'Send method', 'Run', 'Stop', 'Abort', and 'Snapshot'. Below these buttons, there is a field for 'Chromatogram File Name (std 3\_22-III-2006\_010)' with a file name pattern '%q\_%D\_%3n' and a 'Data Recovery...' button. At the bottom, there are 'OK', 'Cancel', and 'Help' buttons.

**Fig. 11. Single Analysis dialog**



## 4.7 Sequence window



The window for defining sequential measurement of several samples is displayed using the  icon in the **Instrument** window.

Measuring is performed by rows in the table. Each row can define measuring of a greater number of samples by the same method or multiple measurements of one sample.

Agilent 1100 - Sequence Ethanol in blood

	Sts.	Run	SV	EV	IV	Sample ID	Sample	Sample Amount	ISTD Amount	Sample Dilut.	Inj. Vol. [µl]	File Name	Std	Lvl	Method Name	Report Style	Open	Open Calib.	Print
1			1	1	1	blank		0,000	0,000	1,000	2,000	Blank	No		Ethanol in blood		<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
2			2	2	1	std1	0.4	0,000	0,000	1,000	2,000	Cal-04	Yes	1	Ethanol in blood		<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
3			3	3	1	std2	0.8	0,000	0,000	1,000	2,000	Cal-09	Yes	2	Ethanol in blood		<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
4			4	4	1	std3	1.4	0,000	0,000	1,000	2,000	Cal-14	Yes	3	Ethanol in blood		<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
5			5	5	1	std4	1.9	0,000	0,000	1,000	2,000	Cal-19	Yes	4	Ethanol in blood		<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
6			6	6	1	std5	2.4	0,000	0,000	1,000	2,000	Cal-24	Yes	5	Ethanol in blood		<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
7			6	6	1	std5	2.6	0,000	0,000	1,000	2,000	Cal-26	Yes	6	Ethanol in blood	Calibration	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
8			7	7	1	0442		0,000	0,000	1,000	2,000	Pers01A	No		Ethanol in blood	Chromatog...	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
9			8	8	1	0445		0,000	0,000	1,000	2,000	Pers01B	No		Ethanol in blood	Chromatog...	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
10																	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>

For Help, press F1      0,24 min. - Running      Vial: 2 / Inj.: 1      File Name: Cal-04      Actix

**Fig. 12. Sequence window**

**ACTIVE** or **PASSIVE** operation is possible. This means that the start-up and duration of individual analyses can be controlled by Clarity station (**ACTIVE** Sequence) or by the auto-sampler and chromatograph (**PASSIVE** Sequence).

### Active sequence:

- Increases reliability of mutual synchronisation between **Clarity**, autosampler and chromatograph.
- Is essential for controlling selected autosamplers using **AS Control** (p/n A26)
- Requires interconnection of a control signal from the **Clarity** station to the autosampler.

**TIP:** Results of sequential measuring can be displayed in the **Chromatogram** window on the **Summary** tab using the **File - Open Chromatogram from Sequence** command (see ch. 4.4.2 - Summary results table on pg. 24).



## 4.8 Program exit

Shut down the program using the **System - Exit** command from the main **Clarity** window.

Before exiting you will always be informed of all running measurements and unsaved files.

## 4.9 Data update and consistency

The **Clarity** station automatically keeps all data updated and consistent. For example, modification of the baseline is automatically manifested in the **Result table**. A change of any parameter updates all related calculations and displayed results. Therefore, the station does not need to contain such commands as **Integrate**, **Recalculate**, etc.



## 5 First "Analysis"

Description of an exemplary analysis will demonstrate basic practical work procedures with the **Clarity** station. Naturally, this description does not exhaust all of the station's capabilities.


**The First Analysis will guide you through the following steps:**

- Program Start
- Project assignement
- Method development (data acquisition parameters, integration parameters, etc.)
- Data acquisition (Sample description, Start, online data monitoring, etc.)
- Chromatogram modification (Integration, etc.)
- Calibration (create and assign calibration file)
- Creating new Template method
- Reporting (Print)

### 5.1 Program start


#### Comments

One of the **Clarity** station's key advantages is immediate measurement without any lengthy preparation and setting.

You can immediately run an analysis measurement by clicking on the  icon in the **Instrument** window and using pushbutton **Run** in the opened **Single Analysis** dialog.

**Note:** *Measurement can also be started immediately after opening the **Instrument** window by the start signal directly from the chromatograph.*

#### Proceed

- Click the  icon on the desktop or select **Start – Programs – Clarity – Clarity Chromatography Station**.  
**Clarity** window will appear.



- Click on the picture of the gas or liquid chromatograph.

## 5.2 Create new project

### Comments

The default projects WORK1 - 4 are predefined for each Instrument. It is recommended to create your own new project to store your data. The Project serves as a base unit to organize corresponding methods, sequences, chromatograms and calibration files.

**Note:** *There are several DEMO projects that contain example data for selected separation techniques. These can be used for learning the station operation. See the DEMO manual (available from [www.dataapex.com](http://www.dataapex.com) or on the installation CD) for the step by step procedure.*


### Proceed

- Open the **Project Setup** dialog using the **File – Projects** command in the **Instrument** window.
- Click the **New** ① button to open **New Project** dialog.
- Set the name of your new project ② and click **OK**.

**Note:** *The **Project Setup** dialog now displays the new project.*

- Specify the ANALYSIS ③ and CALIBRATION ④ subfolders (we recommend to leave the default names).
- Click **OK**.

New MYPROJECT.PRJ ⑤ file in PROJECT folder and corresponding MYPROJECT ⑥ folder in the CLARITY root folder will be created.



**Fig. 13. Creating new project**

All template method (\*.MET) and sequence (\*.SEQ) files will be stored in the MYPROJECT folder.

Calibration files (\*.CAL) and calibration standards (\*.PRM) will be implicitly stored in CALIB and Chromatograms (\*.PRM) in DATA subfolders.

The COMMON ⑦ folder contains files shared among projects (e.g. report styles).

## 5.3 Method development

### Comments

Method contains parameters of the acquisition, integration, device control, calculations, etc.

The template method opened in the **Instrument** window will be copied to the chromatogram after:

- a) acquiring new chromatogram
- b) reprocessing old chromatograms.

The template method is stored as a \*.MET file in the project directory (by default).

**Note:** *In our example the template methods will be stored in the MYPROJECT folder.*

Changes applied to the chromatogram copy of the method (chromatogram method) affect the respective chromatogram only.



The calibration file (\*.CAL) used for identifying and quantitating the peaks is linked by its name in the template method and consequently in the chromatogram method.

In the next step we will prepare template method for new chromatograms with general acquisition parameters.

Later, after measuring the chromatogram, we will modify further parameters (e.g. integration table) in the chromatogram method.

Finally we will save the changed chromatogram method as a template method.

### 5.3.1 Measurement and evaluation parameters


#### Comments

Certain parameters (voltage range, sampling rate etc.) must be set prior to measurement.

#### Caution!

*Changing of these parameters is blocked during measurement because it would render the measurement void.*

#### Proceed



- Create a new NONAME method using the **File - New Method** in the **Instrument**.
- Place the cursor over the chromatograph icon to change the image to four icons.
- Clicking on one of the circled icons will display the contents of the current template method in the **Method Setup** dialog (see Chapter **4.3 - Method Setup dialog** on pg. 21).
- On the **Measurement** tab, fill in the description of the chromatographic conditions and set up the **Enable Autostop**, **Run Time** and **External Start/Stop** options.
- On the **Acquisition** tab, set up the **Range** and **Sample Rate** according to your detector output and expected peak widths.

#### Note:

*The peaks exceeding the specified range will be cut and the respective area lost, however lower*





range settings will increase the resolution at the baseline.

**Note:** The 10 Hz (10 points per second) **Sample Rate** is appropriate for most packed column separations in GC or HPLC, for capillary GC or CE with peak widths below 2 sec the 25 or 50 Hz should be used.

- On the **Integration** tab, the **Global Peak Width** and **Global Threshold** can be set according to the expected peaks.

**Note:** Default settings can be used in most cases. The other entries are best modified in the **Chromatogram** window on already measured chromatogram using the graphical tools.

### 5.3.2 Create and assign a calibration

#### Comments

- On the **Calculation** tab, use the **New** button to create the new calibration file. This will attach it to the template method.

You can use the **View** button to see and modify the file, especially the settings in the **Calibration Options** dialog (**Calibration – Options** menu). The calculation type and other settings may be amended as well.


- The template method is now ready to acquire the first chromatogram - a calibration standard preferably.

## 5.4 Data Acquisition

### 5.4.1 Sample description

#### Comments


The **Chromatogram File Name** field determines the chromatogram name.


Here you can use the tool for automatic chromatogram naming , e.g. ordinal number using symbol **%n**, entering the date using symbol **%d**, etc.




It is advisable to include all information that may be used for searching for the measured chromatogram later on.

## Proceed



- Click on the  icon and fill in the header of the measured analysis (**Sample ID**, **Amount**, etc.).

**Fig. 14. Single Analysis dialog**

- Write the name of the chromatogram in the **Chromatogram File Name** field (you can use the wizard  icon to add variables).


**Note:** *If a file of the given name already exists, you will be advised accordingly at the end of the measurement.*


## 5.4.2 Signal monitoring

### Comments


It is possible to see the online signal from detectors to verify that the detectors are connected and to check the drift and noise.

## Proceed



- Click the  icon to display actual signal from detectors and check drift and noise (see **Fig. 15** on pg. 36).



- Finally, close the window using the  icon or the **File - Exit** command.

### 5.4.3 Analysis run

#### Comments

Acquisition can be started externally from the autosampler or by pressing the **Run** button.

While running, the current measurement time and **RUNNING** inscription will be displayed in the **Instrument** window in the Status line under the diagram.





#### Proceed


- Run the analysis using the external start or the **Run** command in the **Single Analysis** dialog.







### 5.4.4 Analysis monitoring


#### Comment



Clicking on the  icon starts the **Data Acquisition** window, which allows the monitoring of the analysis progress.

Mark the area to be magnified by pressing and holding the left mouse button and marking the desired area and releasing the mouse button. The original display can be returned using  icon or by double-clicking the left mouse button in the graph area.

Data acquisition can be controlled using following icons:  **Single Run**,  **Sequence Run**,  **Sequence Resume**,  **Snapshot**,  **Abort** and  **Stop**.




**Fig. 15. Data Acquisition window**

The measured chromatogram can be compared to an already completed chromatogram, e.g. solvent, calibration standard, etc. Two types of chromatograms can be displayed in the background:

- Simply by using the **File – Set Background Chromatogram** command select which chromatogram is to be displayed in the background.
- Or, if you defined an automatic subtracted chromatogram in the **Method Setup - Advanced** dialog, using the **File - Show Subtraction Chromatogram** command will display the chromatogram chosen for subtraction.

**Note:** *Any chromatograms in the background are displayed only during data acquisition.*

#### 5.4.5 Analysis stop

The analysis can be stopped either by using the  **Stop** command in the **Single Analysis** dialog or automatically after the passing of the **Run Time** set next to the **Enable Autostop** checkbox in the **Method Setup – Measurement** dialog.


At this moment the measured data are integrated, evaluated and saved.








## 5.5 Chromatogram evaluation

### 5.5.1 Display chromatogram

#### Comments




The completed analysis is automatically displayed in the **Chromatogram** window. Automatic analysis display can be prevented by switching the  symbol next to the  icon to the  symbol.


The **Chromatogram** window can also be displayed manually by clicking directly on the  icon. The required chromatogram can be selected using the  icon.

You can magnify any section of the chromatogram in the same way as in **Data Acquisition** window.

#### Proceed

**Note:** *If the chromatogram had opened automatically, skip the next two steps.*

- Click the  icon in the **Instrument** window to invoke the **Chromatogram** window.
- Use the **File – Open Chromatogram** command to open chromatogram.
- Click the left mouse button, hold it, select the area and then release the mouse button to zoom at the part of the chromatogram of interest.




- Doubleclick in the chromatogram area will restore the initial zoom.



## 5.5.2 Chromatogram processing

### Comments




After stopping an analysis it is advisable to check and optimize the integration parameters.

The key parameters are **Peak Width** and **Threshold**. Their simplified definition is:

- The **Peak Width** parameter determines that all distinctly narrow peaks, compared to the defined value, will not be integrated.
- The **Threshold** parameter determines the noise threshold, which means that only peaks of at least double the height of the defined value are detected.

**Note:**

*In reality, however, the situation is more complex as it always combines the mutual effects of both parameters and the peak shape.*



**Fig. 16. Method Setup - Integration**

Both parameters can be applied to the whole chromatogram (Global) or to a defined interval (Local).

The  **Global (Local) Peak Width** and  **Global (Local) Threshold** interactive functions in the **Chromatogram** window will help



to suggest the initial values based on the narrowest peak of interest (peak width) and baseline part with no peaks of interest (threshold).

## Proceed

- Switch to the **Integration** tab in the **Chromatogram** window to see the amendments to the chromatogram method integration table.
- Use the **Global Peak Width** and **Global Threshold** interactive commands to optimize the integration of peaks.
- Limit the integration only to areas where peaks of interest are present using the **Integration Interval** command

**Note:** *This will simplify the baseline detection around the dead volume dips or baseline shifts.*

The other integration tools can be used to get the desired integration, when it is not possible to do so by the previous parameters (see next chapter).

### 5.5.3 Chromatogram modifications



If you wish to edit the automatically generated baseline, you can manually modify the chromatogram whenever you like by:

- either interactively using icons from the tool-bars located on the left edge of the **Chromatogram** window (or corresponding commands from the submenu **Chromatogram - Baseline, Peak and Integrations**). Modifications are saved in the **Integration table** on the **Integration** tab.
- or explicitly enter and change parameters directly in the **Integration table** on the **Integration** tab.

**TIP:** *The **Integration** sheet in the **Chromatogram** window is used for modifying the current chromatogram. The **Integration** sheet in the **Method Setup** dialog is used for editing the whole method for newly measured chromatograms.*



### 5.5.4 Saving of a template method


The changes made to the **Integration table** affect only the current chromatogram. To use those parameters for new chromatograms or reprocess, the current chromatogram method has to be saved as a template method using the command **Method - Save as Template**.

**Note:** *In case the used template method was already saved, you may need to save it under different name, as it is not possible to overwrite method currently opened in the **Instrument** window (or elsewhere).*

## 5.6 Calibration

### Comments

The calibration file (\*.CAL) used for identifying and quantitating the peaks is referred by its name in the template and consequently in the chromatogram method.



Creation and use of calibration curves is an extensive part of the **Clarity** station. Calibration curves created from calibration standards are saved in separate calibration files. Each file can contain virtually an unlimited number of curves calibrated in up to twenty concentration levels plus blank.

Acquisition of calibrated results is done by creating calibration curves, saving them in a calibration file and subsequent linking of a calibration file to a chromatogram, and setting the required type of calculation (external or internal standard).



**Note:** *The station also allows for calibration at several levels and multiple recalibration by optional algorithm. Once again, all automatically or by manual selection of individual compounds, including interactive filling in or control of all data. For all compounds you can also select the calibration curve fit, zero type, type of response, size of identification windows and improved identification by method of reference peaks.*





### 5.6.1 Opening a calibration file

When a calibration file was specified in the template method used for acquiring the chromatogram, you can open it from the **Chromatogram** window.

- Use the **View** button in the **Chromatogram - Results** window.
- Use the  **File - Open Standard** (yellow icon) command to open an integrated chromatogram representing a calibration standard.
- Click the **Add All**  (blue line) icon to automatically transfer available data on each peak to the calibration table.
- In the displayed calibration table, modify compound names in the first column ①.

**Note:** *Table fields can be edited after clicking on them with the mouse.*

Used	Compound Name	Reten. Time	Left Window	Right Window	Peak Type	Peak Color	LOD	LOQ	R B	Resp. Factor	Response	Amount	Resp. Fact	Rec No
6	<input checked="" type="checkbox"/> fructose	6,577	0,100	0,100	Ordrrr		0,000	0,000	A	0,0000	0,0000	0,600	0,0000	0
7	<input checked="" type="checkbox"/> succinic	8,177	0,100	0,100	Ordrrr		0,000	0,000	A	0,0000	31,5600	0,239	0,0076	0
8	<input checked="" type="checkbox"/> acetic	8,550	0,100	0,100	Ordrrr		0,000	0,000	A	0,0000	51,3650	0,600	0,0117	0
9	<input type="checkbox"/> glycerol	8,900	0,100	0,060	Ordrrr		0,000	0,000	A	0,0000	0,0000	1,223	0,0000	0
10	<input checked="" type="checkbox"/> acetic	10,337	0,100	0,100	Ordrrr		0,000	0,000	A	0,0000	24,9850	0,202	0,0081	0
11	<input checked="" type="checkbox"/> methanol	12,710	0,100	0,100	Ordrrr		0,000	0,000	A	0,0000	0,0000	0,104	0,0000	0
12	<input checked="" type="checkbox"/> ethanol	14,833	0,100	0,100	Ordrrr		0,000	0,000	A	0,0000	0,0000	10,647	0,0000	0

- Fill in the known quantity for each compound in the **Amount** column ② under the number of the current calibrated level.

**Note:** *If you have performed all steps as described, you should see a straight-line single-point calibration curve after clicking on the tab of any individual compound.*

- Finally, save the created calibration curves using the **File - Save** command.

### 5.6.2 Link the calibration file with chromatogram or template method

#### Comments




The calibration can be attached (applied) to a chromatogram by selecting its file in the **Chromatogram - Results** window (in **Calibration File (Peak Table)** field).




It can be also attached to a template method file to be applied to new chromatograms (or to reprocess existing ones) in the **Method Setup – Calculation** dialog.

### 5.6.3 Calibration curve use


The simplest way of verifying the correctness of the calibration is to link it to the chromatogram from which it was created.

- Move to the **Chromatogram** window (e.g. by clicking on  icon).
- Here, use the  icon to display the dialog for chromatogram selection.
- Open the chromatogram from which you created the calibration curves. As calibration standards are mostly saved in subdirectory CALIB, first select this calibration subdirectory (e.g. the  (red) icon in the top-right corner of the **Open Chromatogram** dialog). Here, select the chromatogram from which you created the calibration curves.



Then set the calculation method:

- Display the **Chromatogram - Results** tab.
- In the right-hand section, use the **Set** button to select your calibration file and select **ESTD** in the **Calculation** field.



Now the results table should show known quantities for compounds filled in by you, including names and percentage ratios.

In the case of peaks for which no calibrated compound was found based on conformity of



the retention time, the quantity will be zero. If the calibration compound is found and allocated, but no calibration curve exists, the **Peak Type** column will display the word *Error*.


### Caution!

*The chromatogram typically uses calibration curves from the calibration file whose name is shown in the right-hand section of **Results** tab of the **Chromatogram** window. Therefore, this is calibration by reference. Its advantage lies in easy recalibration and changes in case of use of a calibrated file for more chromatograms. Nevertheless, all calibration curves are simultaneously saved directly in the chromatogram, including their history. Results according to these saved calibration curves (including then valid chromatogram modifications) can be recalled if you select the relevant method in the dialog for chromatogram selection - **Open Chromatogram** in the **Method** field.*

File Name:	2506MULTI.PRIM	Detectors:	<input checked="" type="checkbox"/> 1 <input checked="" type="checkbox"/> 2
File Type:	Chromatogram files (*.prim)		
Method:	25.4. 2003 15:22:08 Recent (Linked Calibration)		
Analyst:	25.4. 2003 15:22:08 Recent (Linked Calibration)	Version:	
SampleID:	25.4. 2003 14:54:53 175.4. 2003 14:53:16	Range/Rate:	

## 5.7 Printing

### Comments



Page Setup	<input checked="" type="checkbox"/> Print
<input checked="" type="checkbox"/> Lab. Header	
<input checked="" type="checkbox"/> Report Header	

To facilitate the printing, all the important windows contain icons and .

The report style file actually defines the layout and contents of the printed report.

Information is divided into subsections, corresponding to the displayed tabs. On each tab you can select a more detailed arrangement.

To include the entire section in the report, check the **Print** checkbox in the top-left corner of the tab. The selected section is marked with a tick on the tab label.


**TIP:** A tab can also be selected by double-clicking the left mouse button on the tab.




By default there are different report styles predefined according to the dialog from which



you have invoked the printing commands. The content relevant in the current window is printed. You can easily customize the report styles according to your preference.

## Proceed



- Click on icon  in the **Chromatogram** window.
- Using the **Printer** command on the right-hand side of the **Report Setup** dialog select and set the printer.
- Modify current report style by selecting options on the individual tabs. The settings can be preserved (saved into the report style) using the **Save** or **Save As** commands.
- Using the  icon you can display a print overview and continue to print using the  icon.

# 6 Connecting Autosamplers

This chapter describes wiring of the autosamplers. The configuration varies depending on the type of chromatograph (GC or LC), sequence mode (**ACTIVE** or **PASSIVE**) and presence of optional control modules for direct control of chromatography equipment.

**Note:** For detailed information on **Clarity** sequence modes see ch. 4.7.

### The typical configurations are:

- AS + GC set – **ACTIVE** sequence
- AS + LC set – **ACTIVE** sequence
- AS – **PASSIVE** sequence (GC or LC)
- AS – **ACTIVE** sequence + AS Control
- AS – **ACTIVE** sequence + AS Control + digital acquisition

#### 6.1.1 AS + GC set – Active Sequence

With GC combined with Autosampler, the sample cycle is typically governed by the GC. With the commonly used temperature gradient, the necessary cool down of the system takes variable



time. The sampler is thus synchronized with the GC by a signal wire (READY), allowing next injection only after the GC gets to the **READY** state. The Autosampler performs the injection and starts the GC using another signal wire (START). Any autosampler that is used in the Clarity **Active Sequence** without an **AS Control** module must be synchronized by cable with the **Clarity** station as well as with a chromatograph. The **IN<sub>n</sub>** starting cable should be plugged into the synchronization output (INJECTION) of an autosampler or GC. The **OUT<sub>n</sub>R** cable should be connected to the synchronization input between GC and autosampler.


**Note:** *Further information about **Active Sequence** is available in chap. 4.3.2.2-3. of the **Clarity User Guide**.*

All commonly used autosamplers may be divided in two groups:

- Autosamplers started by closing the contacts on the input (**READY**).
- Autosamplers started by opening the contacts on the input (READY).

#### **Variant A - started by closing the contacts**

The first scheme displays the wiring of an autosampler that is preparing to inject after closing the contact on its input.



**Fig. 17. Wiring of the autosampler – variant A**

The injection will start only after both serially connected contacts (Clarity and the GC) are closed.

After an injection, the autosampler will close the INJECTION contact and thus the command to start the temperature gradient program will be given. At the same time, the chromatograph will close the START contact and thus the command to start the Clarity station acquisition will be given.

If the chromatograph does not have a **START OUTPUT** contact then the starting cable **IN<sub>n</sub>** must be connected directly to the **INJECTION** output on the autosampler (this way, in fact, parallel to the **START INPUT** contact of the chromatograph).

To have the contact on the **INT7** board (INT5 or U-PAD) opened in the initial state, it is necessary to set the **Output Initial** item in the **Digital Output Control** dialog (accessible from the **Clarity** main window) to **HIGH**.

**Note:** The start output mapping of the **Clarity** to individual digital outputs of the A/D converter can be set in the bottom-right corner of the **System Configuration** dialog (see **Fig. 2** on pg. 15). A detailed description of this can be found in the **Reference Guide**.




### Variant B - started by opening the contacts

In the second scheme there is autosampler wiring that conversely waits for output contacts to be opened. This requires different connection (marked by circle).

The **OUTPUT** and **READY** contacts are connected parallel to each other and the autosampler will start its operation after both contacts are opened.

To have the contact on the **INT7** board (INT5 or U-PAD) closed in the initial state, it is necessary to set the **Output Initial** item in the **Digital Output Control** dialog (accessible from the **Clarity** main window) to **LOW**.



**Fig. 18. Wiring of the autosampler – variant B**

**Note:** The start output mapping of the **Clarity** to individual digital outputs of the A/D converter can be set in the bottom-right corner of the **System Configuration** dialog (see **Fig. 2** on p. 15). A detailed description can be found in the **Reference Guide**.

#### 6.1.2 AS + LC set – Active Sequence


In LC systems, the autosampler typically governs the timings. The eventual pump gradient or detector programs are set independently.



Any autosampler that is used in the Clarity **Active Sequence** without an **AS Control** module must be synchronized with the **Clarity** station by cables. The **IN<sub>n</sub>** starting cable should be plugged into the synchronization output (INJECTION) of an autosampler and the **OUT<sub>n</sub>R** cable plugged into the synchronization input (READY) of an autosampler.

**Note:** Further information about **Active Sequence** is available in the **Clarity User Guide**.

Autosampler wiring in LC set for Active Sequence without using the AS Control module.



**Fig. 19.** Wiring of an autosampler in an LC set

**Note:** The labels on the input and output contacts may vary depending on the type of the autosampler. When using additional devices (Detectors, LC Pumps, etc.) it is recommended to connect these devices independently to other digital outputs of the A/D board. Each device will then need a dedicated row in the **Event table** to be started or stopped.

**Note:** When the detector or pump start inputs are connected in parallel to the Clarity start input, problems are often encountered due to instrument grounding.


### 6.1.3 AS - Passive Sequence (GC or LC)

The autosampler used in the Clarity **Passive Sequence** (both GC and LC sets) does not





utilize the **OUT<sub>n</sub>R** digital output cable. All timings are governed by the chromatograph; Clarity only performs one analysis for each start signal received. All synchronization only includes external start up of data acquisition in **Clarity** using the **IN<sub>n</sub>** starting cable.



**Fig. 20. Wiring of an autosampler in Passive Sequence**

A **Passive Sequence** has to be used e.g. in the sets with Headspace autosamplers (without AS Control module).

**Caution!**

*It is not recommended to use the Passive Sequence together with the Control modules.*

**Note:**

*Detailed information on using **Passive Sequence** can be found in the **Clarity User Guide**.*

#### 6.1.4 AS – Active Sequence + AS Control + A/D card

When using the optional **AS Control** (p/n A26) module, all communication is performed through a separate data cable (usually a serial cable connected to a **COM** port).


**Caution!**

*Refer to the corresponding **Clarity Control** manual for wiring specific to your control module.*

Following scheme corresponds to the set with directly controlled autosampler without using the digital acquisition, when data are acquired using an A/D converter (INT7, U-PAD).

In such a case do **NOT CONNECT** the digital output cable **OUT<sub>n</sub>R**.

The starting cable **IN<sub>n</sub>** MUST BE CONNECTED.




**Fig. 21.** *Wiring of an autosampler with AS Control module + A/D converter*

**6.1.5 AS – Active Sequence + AS Control + digital acquisition**

When using optional AS Control module in combination with digital acquisition detectors (e.g. the Duratec DDT3200 module), the connection will be following:

Wiring of the autosampler, with AS Control module and DIGITAL ACQUISITION.



**Fig. 22.** *Wiring of the autosampler with AS control module and digital acquisition*

**Caution!**

*Refer to the corresponding **Clarity Control** manual for wiring specific to your control module.*



## 7 Troubleshooting

If you will not find your answers here, use the [www.dataapex.com](http://www.dataapex.com) website where the Support menu will navigate you to frequently asked questions (FAQ), Clarity email conference archive or contact to DataApex helpdesk.

### 7.1 Locate your problem:

When troubles occur, the fastest way to find a solution for it is to search for it in the following index via the **Dialog** (window), in which the problem occurred, **Error Messages** that appear or according to utilized **Hardware**. The name of the window/dialog is visible in its header.

**Note:** *Names of individual **Clarity Instruments** appear in the header instead of the common term “Instrument”.*

#### Dialog

- Clarity ..... 52, 53, 54, 55
- Data Acquisition ..... 55
- Instrument ..... 54
- Method Setup ..... 54
- Sequence ..... 54
- Single Analysis ..... 54
- System Configuration ..... 54, 55

#### Error Messages

- DEMO – Keylock test failed ..... 53
- DEMO – missing HW Keylock ..... 52
- DEMO - Trial Expired ..... 54
- DEMO - WRONG S/N ..... 53
- DEMO (in the window header) ..... 54
- Disabled (in the status line) ..... 54
- Installation did not pass Windows Logo testing ..... 8
- Other Error Messages ..... 56
- Simulated (in Data Acquisition) ..... 55

#### Hardware

- Hardware key ..... 52, 53



**Note:** *You may find other Error Messages and solutions for problems connected to particular hardware in its respective manuals.*

## 7.2 Problems at station start

### 7.2.1 DEMO – missing HW Keylock

The Rockey (or Sentinel) protective key (dongle) must be plugged in either a parallel or USB port and its driver must be properly installed. Several problems may cause this error message:

1. *Your protective key may not be properly installed.* In the **Control Panel** select the **System** icon, access the **Device Manager** tab and look for "Universal Serial Bus Controllers" - "**Feitian Rockey4 USB**" item (or "**Feitian Rockey4**" for printer port hardlock). It can be also directly in the root folder. If it is not there:

Solution: By using the ROCKINST.EXE program in the HW\_DRIVERS\ROCKEY subdirectory of the **Clarity** workstation.

2. *You might have inserted the **USB** token before installing the **Clarity** station.* Normally **Clarity** station installs proper drivers. In this case, **MS Windows** probably automatically installed incorrect drivers.

Solution: In **Windows XP** it is recommended to use the **System Restore** (menu **Start – Help and Support**) to uninstall the incorrect drivers.

In older versions of **MS Windows** it is necessary to uninstall the incorrect drivers and manually install the correct one as described above.



3. Your protective key may not be correctly connected.

Solution:

**a) For a dongle connected to a parallel port:**

- Try to print on the connected printer.
- Ensure that the printing cable is not too long (>3m).
- Disconnect the printer cable and connect the dongle directly to the port of the computer.

**b) For a USB dongle, check the following:**

- A **USB dongle** may not be used in **Win NT**
- See whether the **USB port** is working (e.g. try to connect a different device, etc.)
- See whether the **HW driver** is installed. In this case the **green LED** on the dongle should be lit.

**7.2.2 DEMO – Keylock test failed**

Your protective key is probably damaged.

Solution:

Contact technical support to replace damaged key.

**7.2.3 DEMO - WRONG S/N**

The **User Code** of the workstation does not match with the code in the protective **Hardware Key**.

Solution:

Use the **Help – User Code** command in the **Clarity** window to enter the correct user code to unlock your station.

The 16-digit **User Code** can be found on the **installation CD** or on the attached CD package inside one of the guides.

**Note:**

The **User Code** dialog does not distinguish upper case and lower case.

If necessary, contact the manufacturer or your dealer to request this code. You will need to provide the serial number of the workstation.

**TIP:**

Be careful not to mix up “I” with “1” on keyboard.




## 7.2.4 DEMO - Trial Expired

*File CLARITY.SNO is missing from the CLARITY directory or is empty.* It is missing or empty erroneously or your **Clarity** station just ended its trial period.

Solution: Copy the file CLARITY.SNO to the main (CLARITY) folder from the **installation CD** (it can be found in the DEMO subdirectory). After starting up the workstation, enter the correct **User Code** using the **Help – User Code** command.

## 7.2.5 DEMO (in the window header)




If you only see the **Demo** caption in the header of the main **Clarity** window without any further description, you have installed the **Clarity Demo** version.

Solution: Uninstall the Demo version and install the full version of Clarity Software.

# 7.3 Problems upon collection of data


## 7.3.1 Data Acquisition – non-functional

Grey icon  with heading “Disabled” and non-functional command Monitor – Data Acquisition



Other manifestations of this error are also: **Method Setup - Acquisition** tab missing, **Method – Acquisition** command non-functional, **Run, Stop, Abort** commands non-functional in the windows **Single Analysis** and **Sequence**). Possible causes are:

### a) Detector not allocated to instrument:



Solution: In the left-hand list **Setup Control Modules** select the correct detector connected to the A/D card you are using and drag it to the corresponding instrument on the right.




If your A/D card is not in the left-hand list **Setup Control Modules**, add it using the **Add** button and repeat the previous step.

**Note:** See chapter 2.4 - Clarity configuration on pg. 14 for detailed description.

- b) You have a licence purchased for data collection from a smaller number of instruments:

Solution:




Open the **System Configuration** dialog from the **Clarity** window using the command **System - Configuration** command and check the tab of the corresponding instrument – **Instrument x**. If the symbol of the curve in the tab header is blue, this indicates an Instrument not usable for measuring.

Check your serial number S/N for example using the command **Help - About** from the main window **Clarity**.

- c) You are using Clarity Offline or a DEMO version, which does not enable the measurement of chromatograms.

Solution:



Check whether there is a blue line with the title **OFFLINE** displayed in the main **Clarity** window under the symbols of the instruments, or the title **Demo** in the window header.

- d) There are problems with the A/D converter INT7 or U-PAD.

Solution:

Consult the more detailed troubleshooting guide concerning the A/D converter in its corresponding manual.

### 7.3.2 Data Acquisition - Simulated

The title “Simulated” is displayed in the Data Acquisition window



The corresponding instrument only displays the simulated curve (from the file CHANNX.DTA), or is allocated as what is known as a DEMO driver.

Solution:

Open the **System Configuration** dialog from the main **Clarity** window using the **System - Configuration** command and check the tab of the corresponding instrument – **Instrument x**.



From the list of equipment allocated to the instrument, take the **Detector x from DEMO driver** field and draw the correct detector of the A/D card you are using from the list on the left.

If you only have Demo detectors in the left-hand **Setup Control Modules** list and your A/D card is missing, open the **Available Control Modules** dialog and using the **Add** button add it to the configuration of the station. Then repeat the previous step.

### 7.3.3 Other Error Messages

You can find the description of other Error Messages and possible problems and solutions for them in other manuals. Here is the list of known possible Error Messages with reference to their descriptions:

- Board Malfunction (INT7, INT9)
- Cannot create detector (INT7, INT9)
- Cannot find driver file \\.\CSWINT70 (INT7)
- Cannot find driver file \\.\CSWINT91 (INT9)
- Cannot load device driver (U-PAD)
- Cannot find first board (INT7, INT9)
- Cannot find second board (INT7, INT9)
- Card not found (INT7) – only older stations
- Error Occurred During Setup (INT7, INT9, U-PAD)
- Cannot establish communication with DataApex U-PAD (U-PAD)

**Note:** *Some of these error messages can also appear if using another hardware than it is listed above. The remedy for such error message should be the same for any device installed. Actualized versions of Clarity Hardware manuals can be found on the DataApex website ([www.dataapex.com](http://www.dataapex.com)).*





## 7.4 Hardware key

### 7.4.1 Rockey USB dongle (re)installation

- Install **Clarity** software from the CD ROM.

#### Caution!

*Be sure to have suitable drivers at hand before inserting the USB token into the PC slot for the first time!*



- Connect the **USB** dongle to a **USB** port of the computer.
- Use the INSTDRV.EXE file from C:\CLARITY\HW\_DRIVERS\ROCKEY\ folder to (re)install the Rockey drivers.
- If this did not help, try the following:
- After connecting the token, **Windows** will recognize the new **Plug and Play** device and will attempt to install a driver. The **Found New Hardware Wizard** will appear:
- Select: "Search for a suitable driver for my device."
- Select: "Specify a location" and then select the path to the CD drive from which you have installed the **Clarity** program, then the HW\_DRIVERS\ROCKEY subdirectory. If you have installed the workstation before connecting the device, you will also find an identical subdirectory in the main directory of the workstation. The driver is in a file named ROCKEY4USB.SYS and its information file is ROCKEY4USB.INF.
- The rest of the installation will be carried out automatically. In the **Start** menu at lower left corner of the screen select, find **Settings - Control Panel**. Then, selecting the **System** icon, verify that the **Device Manager** tab has the correctly installed "Universal Serial Bus Controllers" - "Feitian Rockey4 USB" item.



### 7.4.2 Rockey LPT dongle (re)installation



If you encounter problems with printer port Rockey hardware key, then:

- Verify that the key is properly inserted in the printer port.
- In filemanager (e.g. Windows Explorer) navigate to C:\CLARITY\HW\_DRIVERS\ROCKEY\
- Run the INSTDRV.EXE file to invoke the **Setup Wizard**.
- Select **Repair** in the first dialog and click **Next**.
- In the second dialog check the "**Install parallel driver**" checkbox.
- Select the "**Detect the parallel business**" radiobutton and click the **Next** button.
- In the third dialog, click the **Complete** button.
- You will be asked to Restart Windows.

### 7.4.3 Sentinel dongle installation

Untill Clarity 2.4.1 the dongle had to be installed and present in the PC only when using either the **Clarity Offline** version or when using digital acquisition (e.g. via control module) without an A/D converter. Later versions always require dongle.

#### Printer port dongle



The drivers are automatically installed during the installation of **Clarity** station.

Plug in the dongle in the the printer port (between the printer port and any potential printer).

#### USB dongle

**Caution!**

*Do not insert the Sentinel USB token before prompted to do so by the Sentinel Installer.*




- Run the SENTINEL PROTECTION INSTALLER 7.3.2.EXE program found in the \HW\_DRIVERS\SENTINEL directory of your **Clarity** CD.
- Connect the **USB** token to a **USB** port on the computer when prompted to do so.



## 7.5 System Files (systeminfo.txt file)

The C:\CLARITY\SYSTEMINFO.TXT file contains valuable diagnostic information. Its contents can also be displayed in the Clarity **Help – About – System Files** dialog.



**Fig. 23.** *Help – About – System Files*

The file contains following information (these are examples of what the listings could look like. The information in the program section are visible in the **Fig. 23**):

### System:

*Microsoft Windows 2000 version 5.1 Service Pack 2 (Build 2600)*

### Registered Files:

*CswInt7.dll 11.05.2006 2.4.2.146 C:\Clarity\*

### Files:

*CswAs300.dll 11.05.2006 2.4.2.146 C:\Clarity\*

*Clarity.exe 11.05.2006 2.4.2.146 C:\Clarity\*

### Program:

*Clarity version 2.4.2.146*

...

The last section, **Program**, shows information on installed parts of the **Clarity** station. It shows the version of **Clarity** and the date of



the build, serial number of the station, number of instruments allowed, extensions available, purchased control modules, type and serial number of the hardware key and list of A/D converters/detectors connected to the computer and configured in the station.

The registered files entries should match the installed files in version and location. When there are some discrepancies, it may (but need not) cause some problems.

Solution: If problems occur, use the Remove Clarity item from Clarity group in Windows Start menu. After that, reinstall Clarity.

## 7.6 Sleep Mode

Active **Clarity** station (with **Instrument** window opened) prevents PC from entering sleep mode. This is intentional; otherwise, **Clarity** will not be able to ensure reliable data acquisition.

However, certain types of BIOS may cause problems when the PC enters sleep mode even when the Clarity **Instrument** window is closed. In such case it is recommended to disable the Power Saving features in both Windows OS (for all users) and BIOS.

## 7.7 Switching Users in Windows OS

Switching between User profiles in Windows may freeze the system.

This is caused by problems with communication between the A/D card and system kernel. It is recommended not to switch users on computer where **Clarity** is running.